

SINGLE USE ANORECTAL STAPLERS



ATLAS MEDICAL

The Health Group

ANORECTAL INSTRUMENTS

TCYG Series




TONCARE
Medical Instrument
Manufacture Center

CE 0197 / ISO 13485

ADVANCING FOR LIFE

SINGLE USE ANORECTAL STAPLERS



Increased Cartridge volume

Increased effective volume to accommodate more tissue

Indication Window

Large indication window for precise control of firing

Long Anvil Shaft

Enough length for easier purse-string

Firing Feedback

Provide sound and tactile feedback after firing

Product Overview

TCYG Series has application in selective resection of superior hemorrhoidal mucosa

Single Opening



Suturing-aid



Three Openings



Double Openings



Anal Guide



Anal dilator



Product Code	Color	Diameter(mm)	Number of Stapler	Staple Height (mm)
TCYG-32	White	32	30	4.0
TCYG-34	Green	34	32	4.0
TCYG-36	Green	36	34	4.0



ATLAS MEDICAL

The Health Group



Gallagher Place North U-14
609 Richards Road
Halfway House X94
Midrand
1685



+27 (11) 312 0096



www.atlasmedical.co.za



sales@atlasmedical.co.za

SINGLE USE
ANORECTAL
STAPLERS