

CLOSED SUCTION CATHETERS

Our Closed Suction Catheter is a multiple-use design effective in decreasing the infection of lower respiratory tract infection and can reduce the mortality rate of nosocomial pneumonia.

We provides 24 Hours system and 72 Hours system.



Clinical Practice:

- **Hospital:** economic, no need to stop machine, saving hospital's cost.
- **Nurse:** improving mechanical ventilation's effectiveness, improving airway management quality and successful rescue rate, and lessening the workload of doctors and nurses.
- **Patient:** decreasing the possibility of airway infection, lessening complication, improving patient's life quality, preventing relative pneumonia of respiratory machines

FEATURES AND BENEFIT:

- Double Swivel Elbow is design for endotracheal intubated patients, that provides flexibility without connectors.
- One-way irrigation port to prevent spayback for optimum cleaning.
- Conical waterproof joints for cleaning of the external part of the tube.
- Protective polyurethane sleeve impermeable is flexible and resistant for protection of the healthcare.
- Multiperforated catheter with non-traumatising tip prevents injury.
- 2 side holes for optimum aspiration
- Depth marks enables controlled and deep aspiration.
- Black mark indicates the correct positioning of the tube tip.
- Secured thumb control valve for locked aspiration system.
- Available in difference sizes.

- 24 Hours

REF	Size	ID	Length
PA3524E06	6Fr	2.0mm	45cm
PA3524E08	8Fr	2.7mm	45cm
PA3524E10	10Fr	3.3mm	60cm
PA3524E12	12Fr	4.0mm	60cm
PA3524E14	14Fr	4.7mm	60cm
PA3524E16	16Fr	5.3mm	60cm

- 72 Hours

REF	Size	ID	Length
PA3572E06	6Fr	2.0mm	45cm
PA3572E08	8Fr	2.7mm	45cm
PA3572E10	10Fr	3.3mm	60cm
PA3572E12	12Fr	4.0mm	60cm
PA3572E14	14Fr	4.7mm	60cm
PA3572E16	16Fr	5.3mm	60cm

SUCTION CATHETERS

- Have different connectors

1. PVC Suction Catheter

FEATURES AND BENEFITS:

- Made of medical grade PVC, non-toxic.
- Open tip with 2 lateral eyes.
- Available with different length.
- Color coding for size identification.



- Different types of connectors are available:



Type 1



Type 2



Type 3



Type 4



Type 5

Size	Color Code	Type 1 Control Vent Connector	Type 2 Conical Connector	Type 3 Control Vent Connector	Type 4 Control Vent Connector	Type 5 Control Vent Connector
5Fr	Gray	PA070105	PA070205	PA070305	PA070405	PA070505
6Fr	Light Green	PA070106	PA070206	PA070306	PA070406	PA070506
8Fr	Blue	PA070108	PA070208	PA070308	PA070408	PA070508
10Fr	Black	PA070110	PA070210	PA070310	PA070410	PA070510
12Fr	White	PA070112	PA070212	PA070312	PA070412	PA070512
14Fr	Green	PA070114	PA070214	PA070314	PA070414	PA070514
16Fr	Orange	PA070116	PA070216	PA070316	PA070416	PA070516
18Fr	Red	PA070118	PA070218	PA070318	PA070418	PA070518
20Fr	Yellow	PA070120	PA070220	PA070320	PA070420	PA070520