

1. Guedel Airway (Oral Pharyngeal Airway)



FEATURES AND BENEFITS:

- Made of non-toxic PE, for single use only.
- With color-coded.
- Meet the requirements of ISO 5364.
- Suitable for neonate, infant, child and adult.

REF	Size	Length	Color-coded
PA0401 04	000#	40mm	Pink
PA0401 05	00#	50mm	Blue
PA0401 06	0#	60mm	Black
PA0401 07	1#	70mm	White
PA0401 08	2#	80mm	Green
PA0401 09	3#	90mm	Yellow
PA0401 10	4#	100mm	Red
PA0401 11	5#	110mm	Light Blue
PA0401 12	6#	120mm	Orange

2. Berman Airway (Oral Pharyngeal Airway)

FEATURES AND BENEFITS:

- Made of non-toxic PE, for single use only.
- Clear or color-coded.
- Suitable for neonate, infant, child and adult.

REF	Size	Length
PA0402 04	000#	40mm
PA0402 05	00#	50mm
PA0402 06	0#	60mm
PA0402 07	1#	70mm
PA0402 08	2#	80mm
PA0402 09	3#	90mm
PA0402 10	4#	100mm
PA0402 11	5#	110mm



REF	Size	Length	Color-coded
PA0403 04	000#	40mm	Pink
PA0403 05	00#	50mm	Blue
PA0403 06	0#	60mm	Black
PA0403 07	1#	70mm	White
PA0403 08	2#	80mm	Green
PA0403 09	3#	90mm	Yellow
PA0403 10	4#	100mm	Red
PA0403 11	5#	110mm	Orange



AIRWAYS MANAGEMENT

- Guedel Airway, Berman Airway, Naso-Pharyngeal Airway, Guedel Set

3. Guedel airway Complete Set with Storage Case



FEATURES AND BENEFITS:

- Made of non-toxic PE, for single use only.
- With color-coded.
- Meet the requirements of ISO 5364.
- Suitable for neonate, infant, child and adult.

REF	Product Description
PA3301	Include 8 sizes Guedel airway, 40/50/60/70/80/90/100/110mm

4. Naso-Pharyngeal Airway: Blue Trumpet type and Transparent Flange type

FEATURES AND BENEFITS:

- Made of medical soft PVC, non toxic, Latex Free.
- Smooth surface permits easy insertion and withdrawal.
- Rounded tip for increased patient comfort.
- Suitable for neonate, infant, child and adult.

Blue Trumpet type REF	Transparent Flange type REF	Size ID	Size FR
PA1901 30	PA1902 30	3.0mm	14Fr
PA1901 35	PA1902 35	3.5mm	16Fr
PA1901 40	PA1902 40	4.0mm	18Fr
PA1901 45	PA1902 45	4.5mm	20Fr
PA1901 50	PA1902 50	5.0mm	22Fr
PA1901 55	PA1902 55	5.5mm	24Fr
PA1901 60	PA1902 60	6.0mm	26Fr
PA1901 65	PA1902 65	6.5mm	28Fr
PA1901 70	PA1902 70	7.0mm	30Fr
PA1901 75	PA1902 75	7.5mm	32Fr
PA1901 80	PA1902 80	8.0mm	34Fr
PA1901 85	PA1902 85	8.5mm	36Fr
PA1901 90	PA1902 90	9.0mm	38Fr



- Blue Trumpet type



- Transparent Flange type

FORMED ORAL-PHARYNGEAL AIRWAYS / NASO-PHARYNGEAL AIRWAYS

- For quick and safe airway patency without intubation the user can choose from an extensive range of our oro-pharyngeal airways and naso-pharyngeal airways.
- All versions are developed specifically to suit the customer's characteristic requirements and thus represent a safe and adequate solution for most indications.
- FORMED offers a large number of standard versions made from a variety of materials as well as special models for endoscopy with simultaneous fixation of the tracheal tube.

CLOSED SUCTION CATHETERS

Our Closed Suction Catheter is a multiple-use design effective in decreasing the infection of lower respiratory tract infection and can reduce the mortality rate of nosocomial pneumonia.

We provides 24 Hours system and 72 Hours system.



Clinical Practice:

- **Hospital:** economic, no need to stop machine, saving hospital's cost.
- **Nurse:** improving mechanical ventilation's effectiveness, improving airway management quality and successful rescue rate, and lessening the workload of doctors and nurses.
- **Patient:** decreasing the possibility of airway infection, lessening complication, improving patient's life quality, preventing relative pneumonia of respiratory machines

FEATURES AND BENEFIT:

- Double Swivel Elbow is design for endotracheal intubated patients, that provides flexibility without connectors.
- One-way irrigation port to prevent spayback for optimum cleaning.
- Conical waterproof joints for cleaning of the external part of the tube.
- Protective polyurethane sleeve impermeable is flexible and resistant for protection of the healthcare.
- Multiperforated catheter with non-traumatising tip prevents injury.
- 2 side holes for optimum aspiration
- Depth marks enables controlled and deep aspiration.
- Black mark indicates the correct positioning of the tube tip.
- Secured thumb control valve for locked aspiration system.
- Available in difference sizes.

- 24 Hours

REF	Size	ID	Length
PA3524E06	6Fr	2.0mm	45cm
PA3524E08	8Fr	2.7mm	45cm
PA3524E10	10Fr	3.3mm	60cm
PA3524E12	12Fr	4.0mm	60cm
PA3524E14	14Fr	4.7mm	60cm
PA3524E16	16Fr	5.3mm	60cm

- 72 Hours

REF	Size	ID	Length
PA3572E06	6Fr	2.0mm	45cm
PA3572E08	8Fr	2.7mm	45cm
PA3572E10	10Fr	3.3mm	60cm
PA3572E12	12Fr	4.0mm	60cm
PA3572E14	14Fr	4.7mm	60cm
PA3572E16	16Fr	5.3mm	60cm

SUCTION CATHETERS

- Have different connectors

1. PVC Suction Catheter

FEATURES AND BENEFITS:

- Made of medical grade PVC, non-toxic.
- Open p Åwith 2 lateral eyes.
- Available with different length.
- Color coding for size idenficÅaon.Å



- Different types of connectors are available:



Type 1



Type 2



Type 3



Type 4



Type 5

Size	Color Code	Type 1 Control Vent Connector	Type 2 Conical Connector	Type 3 Control Vent Connector	Type 4 Control Vent Connector	Type 5 Control Vent Connector
5Fr	Gray	PA070105	PA070205	PA070305	PA070405	PA070505
6Fr	Light Green	PA070106	PA070206	PA070306	PA070406	PA070506
8Fr	Blue	PA070108	PA070208	PA070308	PA070408	PA070508
10Fr	Black	PA070110	PA070210	PA070310	PA070410	PA070510
12Fr	White	PA070112	PA070212	PA070312	PA070412	PA070512
14Fr	Green	PA070114	PA070214	PA070314	PA070414	PA070514
16Fr	Orange	PA070116	PA070216	PA070316	PA070416	PA070516
18Fr	Red	PA070118	PA070218	PA070318	PA070418	PA070518
20Fr	Yellow	PA070120	PA070220	PA070320	PA070420	PA070520

Suction Tubing with Yankauer Handle

- Suction connecting tube, Yankauer handle -

1. Suction Connecting Tube

FEATURES AND BENEFITS:

- Made of transparent medical PVC material for better visualization.
- Kink-resistant tooth-like structure to prevent possible occlusion.
- Available with different connectors, which can connect with suction catheter or yankauer handle.
- Length can be customized.

Tube Length	Size: 24Fr, 3/16"	Size: 30Fr, 1/4"
1.8m	PA080118	PA080218
2m	PA080120	PA080220
2.5m	PA080125	PA080225
3m	PA080130	PA080230
3.5m	PA080135	PA080235



2. Yankauer Handle

FEATURES AND BENEFITS:

- Made of transparent material for good visualization during suction.
- The forming angle of handle is designed to make it more convenient and humanized in clinic use.
- Have three types tip: Crown tip, Round tip, Plain tip.

Control Vent	Crown pÅ	Round pÅ	Plain pÅ
with control vent	PA1401	PA1403	PA1405
without control vent	PA1402	PA1404	PA1406

3. Suction Tubing with Yankauer Handle

FEATURES AND BENEFITS:

- Made of transparent medical material.
- Can be combined with different types of yankauer handle.
- Length can be customized.

Tube Length	Size: 24Fr, 3/16"	Size: 30Fr, 1/4"	DescriponÅ
1.8m	PA090118	PA090218	with Yankauer Handle
2m	PA090120	PA090220	with Yankauer Handle
2.5m	PA090125	PA090225	with Yankauer Handle
3m	PA090130	PA090230	with Yankauer Handle
3.5m	PA090135	PA090235	with Yankauer Handle

