



**ATLAS MEDICAL**

The Health Group

## Tracheal and Bronchial Tubes

- Endotracheal, Tracheostomy, Endobronchial -

*"for medical technology, for life care"*

# ENDOTRACHEAL TUBE

- Standard and Reinforced, uncuffed and cuffed

## FEATURES AND BENEFITS:

- Made of medical grade PVC, for single use only.
- Four types: Standard and Reinforced, uncuffed and cuffed.
- Murphy eye reducing the risk of occlusion and maintaining airflow.
- High volume low pressure cuff providing good sealing.
- Radiopaque line allowing clear identification under X-ray.


## 1. Standard Endotracheal Tube

| Size ID | Standard uncuffed | Standard cuffed |
|---------|-------------------|-----------------|
| 2.0mm   | PA030120          | PA030220        |
| 2.5mm   | PA030125          | PA030225        |
| 3.0mm   | PA030130          | PA030230        |
| 3.5mm   | PA030135          | PA030235        |
| 4.0mm   | PA030140          | PA030240        |
| 4.5mm   | PA030145          | PA030245        |
| 5.0mm   | PA030150          | PA030250        |
| 5.5mm   | PA030155          | PA030255        |
| 6.0mm   | PA030160          | PA030260        |
| 6.5mm   | PA030165          | PA030265        |
| 7.0mm   | PA030170          | PA030270        |
| 7.5mm   | PA030175          | PA030275        |
| 8.0mm   | PA030180          | PA030280        |
| 8.5mm   | PA030185          | PA030285        |
| 9.0mm   | PA030190          | PA030290        |
| 9.5mm   | PA030195          | PA030295        |
| 10.0mm  | PA0301100         | PA0302100       |



## 2. Reinforced Endotracheal Tube

| Size ID | Reinforced uncuffed | Reinforced cuffed |
|---------|---------------------|-------------------|
| 3.0mm   | PA030330            | PA030430          |
| 3.5mm   | PA030335            | PA030435          |
| 4.0mm   | PA030340            | PA030440          |
| 4.5mm   | PA030345            | PA030445          |
| 5.0mm   | PA030350            | PA030450          |
| 5.5mm   | PA030355            | PA030455          |
| 6.0mm   | PA030360            | PA030460          |
| 6.5mm   | PA030365            | PA030465          |
| 7.0mm   | PA030370            | PA030470          |
| 7.5mm   | PA030375            | PA030475          |
| 8.0mm   | PA030380            | PA030480          |
| 8.5mm   | PA030385            | PA030485          |
| 9.0mm   | PA030390            | PA030490          |
| 9.5mm   | PA030395            | PA030495          |
| 10.0mm  | PA0303100           | PA0304100         |



# ENDOTRACHEAL TUBE

- Oral Preformed, Nasal Preformed

## FEATURES AND BENEFITS:

- Made of medical grade PVC, for single use only.
- Four types: Oral Preformed and Nasal Preformed, uncuffed and cuffed.
- Murphy eye reducing the risk of occlusion and maintaining airflow.
- High volume low pressure cuff providing good sealing.
- Radiopaque line allowing clear identification under X-ray.

## 3. Oral Preformed Endotracheal Tube

| Size ID | Oral Preformed uncuffed | Oral Preformed cuffed |
|---------|-------------------------|-----------------------|
| 3.0mm   | PA030530                | PA030630              |
| 3.5mm   | PA030535                | PA030635              |
| 4.0mm   | PA030540                | PA030640              |
| 4.5mm   | PA030545                | PA030645              |
| 5.0mm   | PA030550                | PA030650              |
| 5.5mm   | PA030555                | PA030655              |
| 6.0mm   | PA030560                | PA030660              |
| 6.5mm   | PA030565                | PA030665              |
| 7.0mm   | PA030570                | PA030670              |
| 7.5mm   | PA030575                | PA030675              |
| 8.0mm   | PA030580                | PA030680              |
| 8.5mm   | PA030585                | PA030685              |
| 9.0mm   | PA030590                | PA030690              |
| 9.5mm   | PA030595                | PA030695              |



## 4. Nasal Preformed Endotracheal Tube

| Size ID | Nasal Preformed uncuffed | Nasal Preformed cuffed |
|---------|--------------------------|------------------------|
| 3.0mm   | PA030730                 | PA030830               |
| 3.5mm   | PA030735                 | PA030835               |
| 4.0mm   | PA030740                 | PA030840               |
| 4.5mm   | PA030745                 | PA030845               |
| 5.0mm   | PA030750                 | PA030850               |
| 5.5mm   | PA030755                 | PA030855               |
| 6.0mm   | PA030760                 | PA030860               |
| 6.5mm   | PA030765                 | PA030865               |
| 7.0mm   | PA030770                 | PA030870               |
| 7.5mm   | PA030775                 | PA030875               |
| 8.0mm   | PA030780                 | PA030880               |
| 8.5mm   | PA030785                 | PA030885               |
| 9.0mm   | PA030790                 | PA030890               |
| 9.5mm   | PA030795                 | PA030895               |



# ENDOTRACHEAL TUBE

- with Suction Catheter

**Endotracheal Tube with Suction Catheter is designed to provide protection against risk of aspiration and decrease the rate of ventilation-associated pneumonia (VAP) in the hospital.**

## FEATURES AND BENEFITS:

- Made of medical grade PVC, for single use only.
- Two types: Standard cuffed and Reinforced cuffed.
- Murphy eye reducing the risk of occlusion and maintaining airflow.
- High volume low pressure cuff providing good sealing.
- Radiopaque line allowing clear identification under X-ray.

## 5. Endotracheal Tube with Suction Catheter, Standard cuffed

| Size ID | Standard cuffed |
|---------|-----------------|
| 5.0mm   | PA230150        |
| 5.5mm   | PA230155        |
| 6.0mm   | PA230160        |
| 6.5mm   | PA230165        |
| 7.0mm   | PA230170        |
| 7.5mm   | PA230175        |
| 8.0mm   | PA230180        |
| 8.5mm   | PA230185        |
| 9.0mm   | PA230190        |



## 6. Endotracheal Tube with Suction Catheter, Reinforced cuffed

| Size ID | Reinforced cuffed |
|---------|-------------------|
| 5.0mm   | PA230250          |
| 5.5mm   | PA230255          |
| 6.0mm   | PA230260          |
| 6.5mm   | PA230265          |
| 7.0mm   | PA230270          |
| 7.5mm   | PA230275          |
| 8.0mm   | PA230280          |
| 8.5mm   | PA230285          |
| 9.0mm   | PA230290          |



# ACCESSORIES

- Stylet, Holder

## 7. Intubating Stylet

### FEATURES AND BENEFITS:

- Flexible use for assistance of intubation.
- Made of malleable aluminum covered with plastic sheath.
- Have smooth formed tip which reduces trauma of the tissue.
- Optional Accessories: Stylet Adaptor.

#### Intubating Stylet

| REF Number | Size | OD    | Recommended Tracheal tube size |
|------------|------|-------|--------------------------------|
| PA180106   | 6Fr  | 2.0mm | ID 2.5 - 3.5                   |
| PA180110   | 10Fr | 3.3mm | ID 4.0 - 5.0                   |
| PA180112   | 12Fr | 4.0mm | ID 5.5 - 7.0                   |
| PA180114   | 14Fr | 4.7mm | ID 5.5 - 10.0                  |



#### Intubating Stylet, with Adaptor

| REF Number | Size | OD    | Recommended Tracheal tube size |
|------------|------|-------|--------------------------------|
| PA180210   | 10Fr | 3.3mm | ID 4.0 - 5.0                   |
| PA180214   | 14Fr | 4.7mm | ID 5.5 - 10.0                  |



## 8. Endotracheal Tube Holder

### FEATURES AND BENEFITS:

- Provide a fast and secure fixation of the Endotracheal tube.
- Accommodate multiple size of Endotracheal tube.
- Bite block protects patient teeth and gums.
- Allow in-use suctioning of the oropharynx.

| REF Number | Product description      |
|------------|--------------------------|
| PA3601     | Endotracheal tube holder |



# TRACHEOSTOMY TUBE

- Standard and Reinforced, uncuffed and cuffed

**FORMED Tracheostomy Tube conform to the patient's anatomy without pressure or tension.  
The optimized ratio between inner and outer diameter allows the greatest possible inner lumen.**

## FEATURES AND BENEFITS:

- Made of medical grade PVC, for single use only.
- Four types: Standard and Reinforced, uncuffed and cuffed.
- Anatomically shaped 90° design for easy insertion and fixation.
- High volume low pressure cuff providing good sealing.
- Radiopaque line allowing clear identification under X-ray.
- Fixation flange design with simple and safe locking mechanism.
- With inner cannula and comfort neck band.

## 8. Tracheostomy Tube, Standard uncuffed and cuffed

| Size ID | Standard uncuffed | Standard cuffed |
|---------|-------------------|-----------------|
| 3.5mm   | PA220135          | /               |
| 4.0mm   | PA220140          | /               |
| 4.5mm   | PA220145          | /               |
| 5.0mm   | PA220150          | PA220250        |
| 5.5mm   | PA220155          | PA220255        |
| 6.0mm   | PA220160          | PA220260        |
| 6.5mm   | PA220165          | PA220265        |
| 7.0mm   | PA220170          | PA220270        |
| 7.5mm   | PA220175          | PA220275        |
| 8.0mm   | PA220180          | PA220280        |
| 8.5mm   | PA220185          | PA220285        |
| 9.0mm   | PA220190          | PA220290        |



## 9. Tracheostomy Tube, Reinforced uncuffed and cuffed

| Size ID | Reinforced uncuffed | Reinforced cuffed |
|---------|---------------------|-------------------|
| 3.5mm   | PA220335            | /                 |
| 4.0mm   | PA220340            | /                 |
| 4.5mm   | PA220345            | /                 |
| 5.0mm   | PA220350            | PA220450          |
| 5.5mm   | PA220355            | PA220455          |
| 6.0mm   | PA220360            | PA220460          |
| 6.5mm   | PA220365            | PA220465          |
| 7.0mm   | PA220370            | PA220470          |
| 7.5mm   | PA220375            | PA220475          |
| 8.0mm   | PA220380            | PA220480          |
| 8.5mm   | PA220385            | PA220485          |
| 9.0mm   | PA220390            | PA220490          |



# ENDOBRONCHIAL TUBE

- Double-lumen

## 10. Double-lumen Endobronchial Tube


### FEATURES AND BENEFITS:

- Made of medical grade PVC, for single use only.
- Two types: left-sided and right-sided.
- The Pilot balloons printed with Tracheal (Clear) and Bronchial (Blue) for easy identification.
- Smooth inner lumens allow easy passage of Bronchoscopes.
- Pre-loaded stylet to help tube maintain shape for easier intubation.
- Packed with switch connector and two low friction suction catheters.
- Radiopaque line allowing clear identification under X-ray.




\* Packed with Stylet, Switch connector and two Suction catheters.

| Type                  | Size | REF Number |
|-----------------------|------|------------|
| Le <sup>Ä</sup> Sided | 28FR | PA34L28    |
|                       | 32FR | PA34L32    |
|                       | 35FR | PA34L35    |
|                       | 37FR | PA34L37    |
|                       | 39FR | PA34L39    |
|                       | 41FR | PA34L41    |



| Type        | Size | REF Number |
|-------------|------|------------|
| Right Sided | 28FR | PA34R28    |
|             | 32FR | PA34R32    |
|             | 35FR | PA34R35    |
|             | 37FR | PA34R37    |
|             | 39FR | PA34R39    |
|             | 41FR | PA34R41    |



**ATLAS MEDICAL**  
The Health Group

| South Africa | [www.atlasmedical.co.za](http://www.atlasmedical.co.za) | [sales@atlasmedical.co.za](mailto:sales@atlasmedical.co.za) |

