

ENDOSCOPIC LINEAR CUTTING
STAPLER AND RELOAD



ATLAS MEDICAL

The Health Group

ENDO STAPLER

TCNP Series



TONCARE
Medical Instrument
Manufacture Center

CE 0197 / ISO 13485


ADVANCING FOR LIFE

ENDOSCOPIC LINEAR CUTTING STAPLER AND RELOAD




Integrated Machining Anvil

To ensure firm tissue stapling. The thin anvil can also be a non-invasive grasping forcep.



Gecko-Claw Surface Design

To ensure that staples can enter into the groove safely and accurately during clamping and firing.



Droplet Groove Design

Oriented and contact area increased by 30% to ensure better staple formation.



Curved Tip Anvil

Provides visibility and maneuverability around target tissues and vessels.

Product Overview

TCNP Series are intended for resection, transection and anastomosis in abdominal, gynecologic, pediatric and thoracic surgery



Product Code	Shaft Length	Qty/Box
TCNP-S	60mm	3
TCNP-M	160mm	3
TCNP-L	260mm	3

Product Picture	Product Code/Color		Stapling Length	Open Staple Height	Closed Staple Height	Qty/Box
	TCZNX-30B	/White	30mm	2.5mm	1.0mm	6
	TCZNX-45B		45mm			
	TCZNX-60B		60mm			
	TCZNX-30C	/Blue	30mm	3.5mm	1.6mm	6
	TCZNX-45C		45mm			
	TCZNX-60C		60mm			
	TCZNX-30E	/Gold	30mm	4.1mm	1.8mm	6
	TCZNX-45E		45mm			
	TCZNX-60E		60mm			
	TCZNX-30F	/Green	30mm	4.8mm	2.1mm	6
	TCZNX-45F		45mm			
	TCZNX-60F		60mm			
	TCZNX-30H	/Tan	30mm	2.0mm	0.7-1.1mm	6
	TCZNX-45H		45mm			
	TCZNX-60H		60mm			
	TCZNX-30I	/Purple	30mm	3.0mm	1.2-1.6mm	6
	TCZNX-45I		45mm			
	TCZNX-60I		60mm			
	TCZNX-30J	/Black	30mm	4.0mm	1.6-2.1mm	6
	TCZNX-45J		45mm			
	TCZNX-60J		60mm			
	TCZNXK-30B	/White	30mm	2.5mm	1.0mm	6
	TCZNXK-45B		45mm			
	TCZNXK-60B		60mm			
	TCZNXK-30H	/Tan	30mm	2.0mm	0.7-1.1mm	6
	TCZNXK-45H		45mm			
	TCZNXK-60H		60mm			



ATLAS MEDICAL

The Health Group



Gallagher Place North U-14
609 Richards Road
Halfway House X94
Midrand
1685



+27 (11) 312 0096



www.atlasmedical.co.za



sales@atlasmedical.co.za

ENDOSCOPIC LINEAR
CUTTING STAPLER
AND RELOAD