

OXYGEN & AEROSOL THERAPY

Oxygen Mask, Nebulizer kit , Nasal Cannula, Venturi Mask, Silicone Mask
Oxygen Mask with Reservoir Bag / Non Rebreathing Mask

5. Venturi Mask with 6 Diluters



- FEATURES:**
- Made of medical grade PVC, Clear or Green.
 - With 6 diluters for different oxygen concentration: Blue 24%, Yellow 28%, White 31%, Green 35%, Pink 40%, Orange 50%.
 - Locking ring secures flow setting.
 - Includes adaptor for high humidity entrainment.
 - Adjustable nose clip assures comfortable fit.
 - Complete with 7-ft. oxygen supply tubing.

REF	Product Description
PA1301	Pediatric Standard (S)
PA1302	Pediatric Elongated (M)
PA1303	Adult Standard (L)
PA1304	Adult Elongated (XL)

6. Nasal Cannula



- FEATURES:**
- Made of medical grade PVC, Clear or Green.
 - Adjustable elastic band for comfortable.
 - Complete with 7-ft. oxygen supply tubing.

REF	Product Description
PA160101	Standard prongs, Infant
PA160102	Standard prongs, Pediatric
PA160103	Standard prongs, Adult
PA160201	Soft prongs, Neonate
PA160202	Soft prongs, Infant
PA160203	Soft prongs, Pediatric
PA160204	Soft prongs, Adult

7. Nasal Cannula with CO₂/O₂ Sampling tube



- FEATURES:**
- Made of medical grade PVC, Clear or Green.
 - Adjustable elastic band for comfortable.
 - Complete with 7-ft. oxygen supply tubing.

REF	Product Description
PA160301	Standard prongs, Pediatric
PA160302	Standard prongs, Adult
PA160401	Soft prongs, Pediatric
PA160402	Soft prongs, Adult

8. Oxygen Mask with CO₂/O₂ Sampling tube



- FEATURES:**
- Made of medical grade PVC, Clear or Green.
 - Adjustable elastic band for comfortable, snug fit below the ears.
 - Complete with 7-ft. oxygen supply tubing.

REF	Product Description
PA052001	Adult Elongated, XL
PA052002	Adult Standard, L
PA052003	Pediatric Elongated, M
PA052004	Pediatric Standard, S

