

SINGLE USE
EXTRACT BAGS



ATLAS MEDICAL

The Health Group

ENDOSCOPIC ACCESSORIES

TCQW Series



TONCARE
Medical Instrument
Manufacture Center

CE 0197 / ISO 13485

ADVANCING FOR LIFE

SINGLE USE EXTRACT BAGS



Product Overview

TCQW Series has application in collect tissue samples or foreign bodies under endoscopic surgery and removing from the body




High Strength Bag

High strength TPU material with good biocompatibility is used, which is soft and has high strength



Better Support Performance

Provide high quality nickel titanium alloy material for better support performance and fully expansion



Small Diameter Shaft

Pass through the 10 mm trocar, which is convenient for operation



Larger Volume

700ml volume bag is provided to facilitate the retrieval of large specimen at one time

Product Code	Diameter(mm)	L1(mm)	L2(mm)
TCQWA-40	40	/	/
TCQWA-60	60	/	/
TCQWA-80	/	92	82
TCQWA-100	/	113	99
TCQWA-130	/	154	118



ATLAS MEDICAL

The Health Group



Gallagher Place North U-14
609 Richards Road
Halfway House X94
Midrand
1685



+27 (11) 312 0096



www.atlasmedical.co.za



sales@atlasmedical.co.za

SINGLE USE
EXTRACT BAGS