


LAPAROSCOPIC INSTRUMENTS



ATLAS MEDICAL

The Health Group [®]





Monopolar cable connector in the rear of the handle

360° Rotation wheel allows easy access to the tissue

with ratchet control

Ergonomic handle increases comfort during operation




Duckbill




Grasper



Clinch



Rat tooth



Dissector



Babcock



Johan



Curved Scissors



ATLAS MEDICAL
The Health Group



Gallagher Place North U-14
609 Richards Road
Halfway House X94
Midrand
1685



+27 (11) 312 0096



www.atlasmedical.co.za



sales@atlasmedical.co.za